

# AUTOMATIC EXTERNAL DEFIBRILLATOR EOS 1000



Address: 7950 NW R53rd street,  
suite 337, Miami Florida 33166



[www.welchmedical.com](http://www.welchmedical.com)



[ventas@welchmedical.com](mailto:ventas@welchmedical.com)

# Automated External Defibrillators

can be life-saving when cardiac arrest occurs



Remote management



7" HD screen



Portable & lightweight (2.3 kg)



Real-time battery status




Visual and voice instructions guide the rescuer through the entire defibrillation process.




Intelligent biphasic truncated exponential biphasic shock waveform




Turn on the AED



Apply AED patches



Follow voice instructions and deliver a shock



## Features

- High capacity, long life lithium manganese dioxide battery, minimum 200 shocks, 5 years life
- Displayed battery and system status
- Detailed voice instructions and illustrations make rescue easy
- Automated tracking and direct data transfer via integrated 4G module (optional)
- Storage of 24-hour ECG data and up to 1000 events
- Deliver 200 J shock to adults and 50 J shock to pediatrics
- Daily self-checks, timely device status reporting, ideal for remote management

