



# M200 DENTAL CHAIR

Driven by innovation, dedicated to  
perfection



Address: 7950 NW R53rd  
Street, Valley 337, Miami, Florida  
331196



[www.welchmedical.com](http://www.welchmedical.com)



[ventas@welchmedical.com](mailto:ventas@welchmedical.com)

Configuration		
Classification	Configuration	M200
Patient chair	Luxury seat with microfiber leather cushion.	★
	Soft start and stop function	★
	Synchronized chair design	★
	Chair locking system	★
	Spit chair position, treatment chair position, emergency chair position	★
	Anti-collision backrest	★
	Multifunctional motor/pedal	★
Control chair	High-definition LCD screen, with digital display of air pressure and rotation speed.	★
	Integrated LCD negatoscope	★
	3 physicians with 9 sets of chair position memory of the chair.	★
	Intelligent self-test on power-up, error code display	★
Operating tray	Instrument tray water heating function (Can supply hot water to the handpiece and 3-way syringe)	★
	6-position adjustable instrument holder	★
	Handpiece anti-retraction function	★
	Water-air-electricity with an on/off switch control	★
	Independent handpiece water adjustment function	★
	Stainless steel rotating tray	★
	Original anti-fouling silicone pad	★
Main Side Box	Rotating side case	★
	Rotating Ceramic Cuspidor	★
	Imported Water and Air Tube	★
	Independent Water Filter, Dual Filtration Function	★
	Two Water Storage Bottles	★
Assistant Control	3-Way Syringe (Hot)	★
	9-Key Control Panel	★
	Adjustable Strong/Weak Suction	★
	Original Anti-Fouling Silicone Pad	★
Dental Light	LED Sensor Light (8 Imported Bulbs, 3 Light Modes)	★
Floor Box	Recessed Floor Box	★
	Digital Voltage Display	★
	Equipped with a wave filter	★
Dental stool	Professional dental stool (MB1)	★
Auxiliary stool	As an auxiliary stool (MB3)	★
Pipe disinfection	One-button automatic disinfection system	★
	Built-in disinfection module	★
	One-button water source change	★
	Water volume test	★

Options	
Selected Imported Leather	√
Recessed floor	√
High and low speed handpiece	√
Built-in ultrasonic scaler and curing light	√
Endoscope system, endoscope holder	√
Air compressor	√

### Color catalog

Selected luxury microfiber leather seat cushion



Selected microfiber leather seat cushion



Weich Medical is a professional dental equipment manufacturer engaged in the research, development, production, and sales of dental products.



# M200

With disinfection devices



## Digital, intelligent, and personalized

MOS intelligent dental chair control system, for easy exit from work



Intelligent self-test on power-up, error code display.



HD LCD screen, with digital display of air pressure and rotation speed.



9 sets of chair position memory, providing customized chair position settings.



Address: 7950 NW R53rd Street,  
Valley 337, Miami, Florida  
331196



[www.welchmedical.com](http://www.welchmedical.com)



[ventas@welchmedical.com](mailto:ventas@welchmedical.com)

# M200

Opening a new era of induction-controlled disinfection



## Built-in disinfection module

- Full pipe coverage.
- Complete disinfection of all pipes and internal strong and weak suction circuits.
- Built-in cover, using the function once the cover is opened.



## One-button automatic disinfection



- No human intervention: the disinfection program runs automatically
- Offers interruption during disinfection, continuous disinfection after a power outage, automatic low water reminder, and easy emergency management.
- One button to easily exit work: Once automatic disinfection is complete, the dental chair will automatically shut off the water and electricity without any operation, ensuring equipment safety.



## Two-bottle design.

- Water can be poured over the top of the cap, eliminating the need for traditional and tedious bottle unpacking.



## TECHNICAL SPECIFICATIONS

- Maximum chair rise height: 860 mm ± 50 mm
- Minimum chair lift height: 440 mm ± 50 mm
- Continuous working time
- Load capacity: 200 kg
- Dimming mode: 8,000 to 20,000 LUX
- Input voltage: 220V AC 50Hz/60Hz
- Air pressure: 0.55 MPa
- Water pressure: 0.2 MPa – 0.4 MPa



Address: 7950 NW R53rd  
Street, Valley 337, Miami, Florida  
331196



[www.welchmedical.com](http://www.welchmedical.com)



[ventas@welchmedical.com](mailto:ventas@welchmedical.com)

# M200 Dental Chair



### Multifunctional Foot Pedal:

Controls the dental chair's height and the handpiece's operating status, etc. Exceptional appearance, sensitive, and durable.



### Extendable and Rotating Tray:

Stainless steel, sterilizable. Can be rotated 90 degrees horizontally, providing more storage space.



Digital Voltage Display with Voltage Stabilizer: Direct display of the digits to protect the dental chair's circuitry.



### Instrument Tray Water Heating Function:

Can supply warm water to the handpiece and 3-way syringe. Effectively reduces pain during treatment and storage.



Safety Anti-Collision Function: Emergency stop upon resistance to prevent incorrect operation.



### Comes equipped with two stools

Dentist stool, microfiber leather, with lumbar support.

Comfortable Nurse stool, smooth movement, easy to work with.

8 Imported Bulbs

3 Light Modes: White light, yellow light, mixed yellow-white light

2 operating modes: sensor/manual



8-20 Lux	High	100,000
Adjustable brightness	High color rendering: Excellent and comfortable lighting effects	hours of life



Address: 7950 NW R53rd Street,  
Valley 337, Miami,  
Florida 331196



[www.welchmedical.com](http://www.welchmedical.com)



[ventas@welchmedical.com](mailto:ventas@welchmedical.com)