

KLINICO 50



MEDICAL SUCTION MACHINE

FEATURES

- Oil-free piston pump.
- Dual operation via foot pedal and automatic function.
- 5-liter capacity (2.5 liters per bottle)
- Equipped with a valve to stop suction once both bottles are full.
- Fast electric suction unit
- All-ABS panel
- Removable bottle cap with safety valve against overflow.
- Easy to transport, low noise, high negative pressure, and high flow rate.
- Can be used in any surgical procedure.
Absorbs phlegm, thick liquids, is easy to clean, and resistant to disinfectant chemicals used in hospitals.
- Technical Specifications
 - Power Supply: AC220V \pm 22V 50Hz \pm 1Hz
 - Maximum Negative Pressure: \geq 0.06MPa
 - Noise: \leq 60 dB(A)
 - Power: \leq 180VA
 - Pumping Rating: \geq 20 L/min
- Reservoir Capacity: 2500mL per bottle, 2 pieces.
- Pressure Range: Greater than 90kPa (Maximum Vacuum)
- Adjustable Vacuum Range: 20kPa to 90kPa



ADDRESS: 7950 NW 53rd Street, Suite 337; Miami,
Florida 33166
Tel. (305) 440-5155
Fax: (866) 480-9591
Email: ventas@welchmedical.com
Website: www.welchmedical.com