



OPERATION TABLE

PHOENIX 200

Address: 7950 NW 53rd Street, Suite 337
Miami, Florida 33166, USA

www.welchmedical.com
ventas@welchmedical.com

TECHNICAL SPECIFICATIONS

Electric Operation Table Phoenix 200

 **Technology**
THAT SAVES LIVES

Technical Parameters

Dimension	2025x520mm
Min Height	750 mm
Max Height	1050 mm
Sliding movement	300 mm
Back Section	Upward: 70° Down: 30°
Lateral Tilt	Left Tilt: 18° Right: 18°
Trendeleburg	20°
Anti-Trendeleburg	25°
Kidney Bridge	120mm
Head Section	Upward: 60° Down: 90°
Leg Section	Up 40° Down: 90° Outreach: 90°

Technical Configuration

Anaesthetic Screen	1 pcs
Shoulder support	1 set
Arm support	1 set
Body support	1 set
Leg support	1 set
Mattress	1 set

Features

WELDING | SAFETY

Stainless steel base, easy to clean, easy to disinfect



PECULIARITY

The head plate can rise 40 degrees and fall 90 degrees



PRODUCT

Leg plate up 40°, down 90°



DESIGN

A backup controller is attached to the column



MATERIAL

Memory sponge mattress, soft and comfortable




PECULIARITY

Operating table remote control is an important tool to control the operating table, through the remote control can easily control the operation of the operating table, left and right tilt, forward and backward



ATTACHMENT


Comes standard with a pair of arm rests



ATTACHMENT

Comes standard with a pair of leg rests





REMOTE CONTROL

Electric Operating Table

