



SZ-10AW

Mobile Oxygen Generator



FEATURES

Mobile Oxygen Generator SZ-10AW

- Durable, high-traffic housing with wheels for easy transport and side tie-downs.
- Oxygen generation method: Pressure Swing Adsorption (PSA).
- LCD screen with timer function.
- Nebulizer outlet on the front of the unit.
- Precise flow rate control via knob and flowmeter.
- Humidifier reservoir.

TECHNICAL SPECIFICATIONS

Mobile Oxygen Generator SZ-10AW

- Approximate maximum flow rate: 10 L/min
- Power consumption: 10 L ≤700 VA
- Oxygen density up to 90% (±3%)
- Noise level: ≤65 dB(A)
- Power supply: 220V, 50Hz

